

HbA1c Control

Components:

HbA1c

Product description

HbA1c Controls are lyophilized control materials based on human whole blood. Designs are carried out according to our customers' specifications.

Product features

- prepared from human whole blood
- excellent commutability
- 1 week stability after opening

Biohazard

The used human blood donations of European origin are tested by CE-marked kits and are non-reactive for HBsAg, Anti-HIV 1/2 and Anti-HCV. In addition HCV and HIV are tested by PCR.

Volumes. 0.5 ml

Storage and Stabilities

Lyophilized: 2 years at +2/8°C
After reconst.: 1 day at +18/25°C
1 week at +2/8°C

(Refer to packaging insert for exceptions)

Order Information

10123 HbA1c Controls